

TABLET FOUR-USAGE TESTER

Features:

- Intellectualized, auto-control, auto-test, auto-diagnose, auto-alarm, auto-display, auto-memorize, auto-reset, auto-retest.
- The circular water current pumped by magnetism pump can heat the system equally, and the bath liquid can achieve an even temperature.
- Parameters can be preset at any time, the preset and real data are displayed alternately in time.
- Two disintegrations run synchronously. Auto stop at timing point.
- The hardness has high accuracy press sensor and hardness data is displayed in Newton.
- It tests tablet hardness continuously. The tablet is loaded and pressed artificially.



Specifications:

Model	PJ-3
Stations of Dissolution	3
Temperature range	(20 ~ 45) °C ± 0.3 °C
Speed range	(20 ~ 200)rpm ± 2 rpm
Automatic time range	(1 ~ 900) min ± 0.5 min
Paddle radial run-out	± 0.5 mm
Basket radial run-out	± 1.0 mm
Basket assemblies	2
Nacelle back and forth frequency	(30 ~ 32) / min
Nacelle back and forth range	(55±2) mm
Cylinder number	1
Drum' s size radius:	Φ286mm depth: 39mm
Tablet fall height	156 mm
Rotation speed	25 rpm ± 1 rpm
Rotation count	100 r ± 1 r
Hardness range	(2 ~ 200) N
Probe procession range	(3 ~ 20) mm
Dimension (mm)	590×300×480
Power	110/220V/50-60Hz,600W